

CORRIGENDUM
STEREOTACTIC CRANIAL FRAME (TENDER NO:SCT/004/IMP/2017-18/GG)

EXISTING	TO BE READ AS
<ul style="list-style-type: none"> • The system should have a coordinate frame g/other coordinates to be specified 	<p>The system should have all coordinate frames</p>
<ul style="list-style-type: none"> • Disposable inserts to be offered 	<ul style="list-style-type: none"> • Disposable pins or inserts should offered whenever such inserts become essential for day today use.
<ul style="list-style-type: none"> • 1 box with 25 x 4 pieces of sterile inserts to be used with insulated fixation posts. The inserts are delivered sterile and pre-packed four and four in small pouches. 	<ul style="list-style-type: none"> • The inserts are delivered sterile and pre-packed four and four in small pouches.
<ul style="list-style-type: none"> • The reusable fixation screw should ensures quick and rigid skull fixation without compromising imaging quality. Reusable fixation screws are to be made of coated hard aluminum. Each set to consists of 20 pairs of reusable screws in various lengths. 	<ul style="list-style-type: none"> • The reusable fixation screw should ensures quick and rigid skull fixation without compromising imaging quality. Reusable fixation screws are preferably made of coated hard aluminum or light and low artifact material. Each set to consists of 20 pairs of reusable screws in various lengths.
<ul style="list-style-type: none"> • offered system to be a rigid and accurate fixation of coordinate frame g to the operation table headrest. Should have the three point fixation allows for high stability even if the patient is shaking. 	<ul style="list-style-type: none"> • Offered system to be a rigid and accurate fixation of the coordinate frames to the operation table headrest. Should have the three point fixation allows for high stability even if the patient is shaking.